

FORENSIC CHEMISTRY

TOXICOLOGY SECTION BLOOD ALCOHOL TRAINING MANUAL

Blood Alcohol Training Manual: Doc # = 010 Originally issued 5-8-2017 Approved by: Forensic Lab Director – Lauren Niskach Date Revised: 07-19-2019

Page 1 of 20



Introduction

This procedure details the training requirements that must be completed before any new employee in the Toxicology Section may be allowed to analyze case specimens without close supervision (i.e., independent casework). Prior to completion of the modules within this manual new employees may work with case specimens as a part of their training, but may do so only if assisted and supervised by a more experienced analyst or supervisor.

Documented completion of each module within this manual, and completion of the safety training, satisfies all training requirements, and will allow the analyst to obtain the required State Certification for analysis of blood alcohol samples.

Following completion and authorization the Forensic Lab Director, any new analyst is allowed full access to all case specimens within the laboratory and may contribute to the completion of a case in any way necessary, including case / batch technical review.

During the training period, if a trainee's performance is unacceptable in any specific section, the Training Coordinator shall notify the trainee and Forensic Lab Director / Quality Manger that performance is/was unacceptable. Notification of the unacceptable performance will be retained in the training binder along with all documentation indicating the trainee 'repeated the section'.

The Forensic Lab Director/Quality Manager may require the trainee to repeat the failed section, and/or augment the training material in the specific section, and/or take other appropriate action as management deems necessary (i.e., disciplinary). Appropriate documentation will be retained in the training binder. Subsequent failures, of either the original or additional modules during the training period will be brought to the attention of higher management by the Forensic Lab Director/Quality Manager.

MODULE 1:

Administrative Procedures and Orientation

Objectives:

Completion of this module will provide the trainee with relevant knowledge in the following administrative areas:

- Administrative policies and procedures for the State of Maine
- Administrative policies and procedures for the Health and Environmental Testing Laboratory: Forensic Chemistry Section

Blood Alcohol Training Manual: Doc # = 010 Originally issued 5-8-2017 Approved by: Forensic Lab Director – Lauren Niskach Date Revised: 07-19-2019

Page **2** of **20**



Required Tasks:

- 1. New employees will be processed through the State Of Maine / Health and Environmental Testing Labs administrative offices prior to reporting to the Forensic Chemistry Laboratory.
- 2. The trainee will be given a tour of the HETL building, and specifically the Forensic Chemistry Lab, and introduced to lab personnel.
- 3. Assigned Training Coordinator will provide the trainee with an explanation of the purpose of the training program, the expected course of events and outcomes, instructions on maintenance of a training notebook / binder, and a copy of the approved training plan. Trainee will also have the opportunity to ask any questions related to training, and what is expected from all parties.
- 4. Trainee and Training Coordinator will meet with the Forensic Lab Director to discuss items outlined in the completion checklist below, and to answer any questions.
- 5. Trainees will be required to complete reading or viewing of the following materials:
 - a. Laboratory Quality Manual
 - b. Laboratory Standard Operating Procedures Manual
 - c. Relevant sectional Procedural Manuals
 - **d.** Forensic Chemistry Safety Program Manual, which includes sections on blood borne pathogens (SharePoint)
 - e. Chemical Hygiene Plan
 - f. American Academy of Forensic Sciences Code of Ethics
 - g. American Board of Forensic Toxicology Code of Ethics
 - **h.** ANAB Guiding Principles of Professional Responsibility
 - i. Current version of ANAB manual under which the laboratory is currently accredited.
- 6. It should be noted that the Trainee may be working on multiple modules at the same time. Some modules, or requirements within a module may take a short time to complete, but other items may require more extensive practice and time for the trainee to master.



Checklist: Orientation and Administration

1. Trainee has completed Ne	ew Employee Processing	
Trainee	Training Coordinator	Date
2. Trainee has completed St	ate of Maine / HETL New Employe	ee Orientation
Trainee	Training Coordinator	Date
familiar with the location meeting locations. Train Universal Precautions and	TL Safety Program Manual. Traine of fire extinguishers, emergency rounce has discussed with the Forensed the expectations for the proper us ninee understands the HETL policy for the proper uses the expectation of the proper use the expectation of the expectations.	ites and exits, and the evacuation ic Lab Director the concept of se of PPE and safety equipment
Trainee	Training Coordinator	Date
*	led an explanation of the training ce of a training notebook, and a copy	
Trainee	Training Coordinator	Date

- **5.** Trainee has met with the Forensic Lab Director and discussed the following issues:
 - Laboratory mission, structure, organization, and capabilities
 - Types of casework commonly encountered
 - ANAB procedures and oversight
 - Aspects of the Toxicology Laboratory Quality Assurance Program
 - Instrument maintenance and documentation
 - Technical and administrative review of batch and case files
 - Preparation, use, and documentation of calibrators/controls
 - Proficiency Testing Program
 - Forms currently in use in the laboratory

Blood Alcohol Training Manual: Doc # = 010 Originally issued 5-8-2017 Approved by: Forensic Lab Director – Lauren Niskach Date Revised: 07-19-2019

Page 4 of 20



- Hours of operation
- Time and attendance, requests for time off
- Annual and sick leave
- Continuing education and travel requests
- Supply ordering and annual budget process
- Policies regarding confidentiality and reporting of results, including documentation of case communications, telephonic results, and inquiries from the press or other media
- Policies regarding testimony
- Request for pre-trial meetings or depositions in a criminal case
- Request to testify in a grand jury proceeding or preliminary hearings
- Providing lab reports to other agencies
- Discussing the results/conclusions of another analyst's work or the work performed in another section
- Discussing details regarding an on-going investigation, and DHHS policy regarding HIPPA
- Disclosure of case details or comments regarding HETL or HETL Toxicology section via current social networking means (such as: forensic websites, forums, Facebook, internet blogs, chat rooms, etc)
- The potential for re-examination of evidence
- Membership dues and attendance at professional meetings

Trainee	Forensic Lab Director	Date



6.	Trainee 1	nas compl	leted read	ing of	the fo	llowing	material	s:
----	-----------	-----------	------------	--------	--------	---------	----------	----

			Completion Date:	Trainee Initial
	Quality Manual			
	Standard Operating Procedures Manu	ıal		
	Blood Alcohol Procedures Manual			
	American Academy of Forensic Scien Code of Ethics	nces		
	American Board of Forensic Toxicolo Code of Ethics	ogy		
	ANAB Guiding Principles of Professional Responsibility			
	Current ANAB Manual and Board Interpretations			
Train	nee Traini	ng Coordinator	Date	
Admi	inistrative Matters and Procedures de	ocumented, co	mpleted, and revie	ewed.
Train	Traini	ng Coordinator	Date	

Blood Alcohol Training Manual: Doc # = 010 Originally issued 5-8-2017 Approved by: Forensic Lab Director – Lauren Niskach Date Revised: 07-19-2019



Module 2:

Handling and Documentation of Physical Evidence

Objectives:

After completion of this module, the trainee will have the requisite knowledge, skills, and abilities involved in:

- Preservation of hard copy chain-of-custody for physical evidence
- Use of the laboratory LIMS system to document service requests and evidence receipt, and return.
- Compliance with Toxicology Laboratory policies regarding security, handling, packaging, labeling, and preservation of evidence to prevent loss, deterioration, or cross-contamination, including:
 - o Security systems in place
 - Verification and documentation of evidence received
 - o Sealing of evidence
 - o Procedures in the event of receiving unsealed evidence
 - o Procedures in the event of receiving an expired blood kit
 - o Marking of evidence
 - Evidence storage
 - o Final disposition of evidence
- Examination documentation to include:
 - o Initials or signature of analyst and toxicology laboratory case number
 - o Procedures for strikethrough and interlineation
 - o Detailed information on condition and description of the evidence
 - o Analyses/examinations performed
 - o Composition and disposition of case files

Required Tasks:

- The trainer will discuss the topics listed in the objectives as they relate to evidence handling, packaging, labeling, and preservation, chain of custody, and examination documentation.
- The trainee will observe utilization of the laboratory LIMS system to document service requests and evidence transfer.
- The trainee will observe the physical receipt, transfer, storage, and return of evidence, and the corresponding documentation of hard copy chain of custody.
- The trainee will observe the marking or labeling of physical evidence.
- When the trainee is comfortable with and understands the submission process, the training coordinator will pretend to be a customer submitting evidence, and the trainee shall receive said evidence, completing the contract for lab services and chain of custody documents. No labels printed, and this case will not be entered into StarLims.

Blood Alcohol Training Manual: Doc # = 010 Originally issued 5-8-2017 Approved by: Forensic Lab Director – Lauren Niskach Date Revised: 07-19-2019

Page **7** of **20**



• The trainee is required to read the relevant sections of the Quality Manual and Evidence Manual dealing with chain of custody, receipt and return of evidence.

Study Questions (Answer in writing):

- Explain the chain of custody system used by the Toxicology section.
- How is access to the building, the laboratory, and toxicology evidence controlled, secured, and documented?
- Define a proper seal.
- Discuss the proper marking/labeling of evidence
- Who has access to the secured refrigerator in the laboratory?
- What information should be included on each page of your case file?
- How are spelling or other errors handled in note-taking?
- What documentation should be included in a case file?
- How are case files stored within the laboratory?
- How is batch and QC data stored within the laboratory?
- What is the final disposition of evidence when individual casework is completed?



Checklist/Evaluation for Handling of Physical Evidence

1.		ate to evidence handling, pa	or and discussed the topics listed in the ckaging, labeling, preservation, chain of
D	ate	Trainee	Training Coordinator
2.	Trainee has observe documenting on chain		nsferring, and returning evidence and
D	ate	Trainee	Training Coordinator
3.	Trainee has observed/s	assisted in preserving and stor	ring evidentiary materials.
D	ate	Trainee	Training Coordinator
4.	Trainee has observed/if warranted.	assisted in marking of eviden	ace, and description of evidence materials
D	ate	Trainee	Training Coordinator
5.	Trainee has provided s	satisfactory answers to oral stu	udy questions.
D	ate	Trainee	Training Coordinator
6.		ted proficiency in use of the I	LIMS system and completed the exercise evidence.
D	ate	Trainee	Training Coordinator

Blood Alcohol Training Manual: Doc # = 010 Originally issued 5-8-2017 Approved by: Forensic Lab Director – Lauren Niskach Date Revised: 07-19-2019

Page **9** of **20**



7. Tr	rainee has completed the require	ed reading of the	following:	
N	Manuals to Read:		Completion Date:	Trainee Initial
E	Evidence Manual			
	Date	Trainee	Training	g Coordinator
Hand	ling and Documentation of Pl	nysical Evidence	e documented, comple	eted, and reviewed.
	Date	Trainee	Trainin	g Coordinator



Module 3:

Blood Alcohol

Objectives:

After completion of this module, the trainee will have the requisite knowledge, skills, and abilities involved in:

- GC Headspace maintenance and preparation for the run
- Sample check-in to validate accurate information such as names, case #, etc
- Prepare Internal Standard Solution used in Blood Alcohol assay
- Prepare Calibration curve and completing associated forms
- Sample, standard and control preparation
- Developing a sequence list and running samples
- Analyzing the data generated and completing associated forms
- Uncertainty of Measurement
- Creating final report for customer

Required Reading:

- Blood alcohol procedures manual.
- Stability of Ethanol in Blood and Urine Samples: Slavka Mandic-Radic, Gordana Dzingalasevic, Nevena Lukovic. Journal of Molecular Biology, 2007: 26 (3)
- Stability of Ethanol in Human Whole Blood Controls: An Interlaboratory Evaluation. Dubowski, Kurt M., Gadsden, Sr., Richard H., Poklis, Alphonse. Journal of Analytical Toxicology, Vol. 21, October 1997.
- Long-Term Blood Alcohol Stability in Forensic Antemortem Whole Blood Samples. Tiscione, Nicholas B., Vacha, Ruth E., Alford, Ilene., Yeatman, Dustin Tate, and Shan, Xiaoqin. Journal of Analytical Toxicology, 2015;39: 419-425.
- The Effect of Temperature on the Formation of Ethanol by *Candida Albicans* in Blood. Chang, Joyce., and Kollman, S. Elliot. Journal of Forensic Sciences, Vol 34, No. 1, Jan. 1989, pp.105-109.
- Comparison Among Plasma, Serum, and Whole Blood Ethanol Concentrations: Impact of Storage Conditions and Collection Tubes. Penetar, David M., McMeil, Jane F., Ryan, Elizabeth, and Lukas, Scott. Journal of Analytical Toxicology, Vol. 32, Sept. 2008.
- Garriott's Medicolegal Aspects of Alcohol, (5th Edition or later). Due to the size of this volume, it is not required that the trainee read every page of every chapter at this time. However, it is important that the trainee review the book and gain a simplistic understanding regarding the value of this one book and that it addresses almost all aspects of alcohol related issues. The trainee is encouraged to spend free time reading this book.

Blood Alcohol Training Manual: Doc # = 010 Originally issued 5-8-2017 Approved by: Forensic Lab Director – Lauren Niskach Date Revised: 07-19-2019

Page **11** of **20**



Required Tasks:

- The trainer will discuss the topics listed in the objectives as they relate to the blood alcohol assay
- The trainee will observe the entire process as it is performed by the trainer (calibration, case preparation and analysis, completion of all associated forms, entering data into LIMS, and report writing.
- The trainee is required to read and follow the Blood Alcohol Procedures and demonstrate the ability to follow that method and generate results.
- The Trainee will prepare a lot of Internal Standard Solution with the trainer and test that it is suitable for use.
- The Trainee with discuss PerkinElmer instrument maintenance requirements and observe how to perform maintenance tasks.
- The Trainee will practice pipetting using the Hamilton Diluter until trainee feels comfortable using the device.
- After watching the calibration procedure, the Trainee will prepare and successfully complete a minimum of 5 Calibration Curves with minimal assistance from the trainer. Trainee will complete all paperwork / forms and submit to Trainer for review.
- The Trainee will prepare and successfully complete a full volatiles calibration curve, complete all paperwork/forms, and submit to Trainer for review.
- The trainee will watch the Trainer complete the process of sample preparation and analysis, including all associated QC samples (standards, whole blood and serum controls), completion of all paperwork and forms, preparation of individual case folders, and writing of reports in StarLims.
- After watching the Trainer complete the analysis process, the Trainee will successfully complete the process of sample preparation and analysis using proficiency test samples reserved for training purposes. Trainee will successfully complete <u>at least 3</u> batches (each batch having a minimum of 36 samples/controls), with all associated QC samples, complete all paperwork/forms, and submit to Trainer for review. The first batch performed by the Trainee will be observed by the Trainer, and the remaining subsequent batches will be performed independently. Mock reports will be created for 3 samples, having the Trainee demonstrate the ability to report the correct result and calculate the UoM for that result.
- The trainee is required to watch and assist in the creation of BAC reports using StarLims. There is no minimum established as it is recognized that mastering the Starlims software package may be quick, or lengthy depending on the computer literacy of the trainee.

Study Questions (Answer in writing):

• Explain the process of GC Headspace maintenance, including frequency.

Blood Alcohol Training Manual: Doc # = 010 Originally issued 5-8-2017 Approved by: Forensic Lab Director – Lauren Niskach Date Revised: 07-19-2019

Page **12** of **20**



- How are samples checked once analysis starts?
- Describe the pipetting process and the use of the Hamilton Diluter
- Describe the process of creating a sequence and running samples.
- Describe the requirements a calibration curve must meet before being used for casework
- Describe the requirements / acceptance criteria for duplicates and replicates of casework samples
- Describe the theory of uncertainty of measurement, and the specific uncertainty of measurement that is applied to blood alcohol samples?
- How is data reviewed and then how are results entered into StarLims



Checklist/Evaluation for Blood Alcohol Analysis

	met with the training coordinate relate to Blood Alcohols.	or and discussed the topics listed in the
Date	Trainee	Training Coordinator
	he Blood Alcohol Procedures Ma location of Garriott's Medicolega	nual and other required readings, along al Aspects of Alcohol within the
Date	Trainee	Training Coordinator
	-	preformed at least 5 calibrations that meet orms/paperwork for a complete calibration
Date	Trainee	Training Coordinator
11. Trainee has comp trainer.	leted the Blood Alcohol process v	while under close supervision of the
Date	Trainee	Training Coordinator
12. Trainee has succe	ssfully prepared and analyzed 3 b	atches.
Date	Trainee	Training Coordinator
13. Trainee has provid	ded satisfactory answers to oral/w	ritten (circle one) study questions.
Date	Trainee	Training Coordinator

Blood Alcohol Training Manual: Doc # = 010 Originally issued 5-8-2017 Approved by: Forensic Lab Director – Lauren Niskach Date Revised: 07-19-2019

Page **14** of **20**



14. Trainee has demonst	trated proficiency in use of the LIMS sy	stem, including the creation of
reports.		
Date	Trainee	Training Coordinator
Blood Alcohol: comple	ted, and reviewed.	
Trainee	Training Coordinator	Date



MODULE 4:

Criminal and Civil Law, Procedures, and Testimony

Objectives:

Completion of this module will provide the trainee with relevant knowledge in the following areas of criminal and civil law and procedures:

- The structure of the court system in the state of Maine
- Current Maine Rules of Evidence concerning relevancy of forensic evidence and expert testimony
- Current Maine statues regarding impaired driving
- Brady and Giglio material, analyst malfeasance

Completion of this module will also provide the trainee with the requisite knowledge, skills, and abilities to provide effective testimony regarding the analysis of forensic evidence, including courtroom procedures involving:

- Oath and subsequent direct examination, including
 - Oualifications and voir dire
 - o Recognition and identification of evidence
 - o Chain of Custody
 - o Toxicology reports
 - o Explanation of examinations and conclusions
- Cross examination and redirect
- Use of visual aids or presentation materials

Required Reading Materials:

The trainee is required to complete reading or viewing of the following materials / websites related to The Criminal Courts within the State of Maine and Brady / Giglio issues.

- Maine Court Systems (by County): http://www.courts.maine.gov/maine courts/
- Current Rules of Evidence in Maine: http://www.courts.maine.gov/rules_adminorders/rules/text/mr_evid_2015-9-1.pdf
- Maine OUI Laws: http://www.maine.gov/dps/bhs/impaired-driving/laws.html
- Brady Disclosure: https://en.wikipedia.org/wiki/Brady_disclosure
- Practical Guide to Brady Motions: http://www.ncids.org/Defender%20Training/2008%20New%20Felony%20Defender%20Training/BradyHandout.pdf
- Bullcoming vs. New Mexico: https://www.supremecourt.gov/opinions/10pdf/09-10876.pdf

Blood Alcohol Training Manual: Doc # = 010 Originally issued 5-8-2017 Approved by: Forensic Lab Director – Lauren Niskach Date Revised: 07-19-2019

Page **16** of **20**



- Mitchell V. Wisconsin: https://www.supremecourt.gov/opinions/18pdf/18-6210_2co3.pdf
- Giglio vs. United States: https://en.wikipedia.org/wiki/Giglio v. United States
- U.S. Justice Department Policy on Giglio: https://www.justice.gov/archives/ag/policy-regarding-disclosure-prosecutors-potential-impeachment-information-concerning-law

Required Tasks:

- 1. Observe a minimum of 2 courtroom testimonies for a court qualified examiner in the discipline, and one testimony of a qualified examiner from another discipline within the Forensic Chemistry Section. Discuss the testimony with the appropriate examiner/analyst.
- 2. Trainee will meet with the Forensic Lab Director / Quality Manager after reviewed the required reading to discuss violations of Brady / Giglio, as well as general analyst malfeasance in the Forensic Sciences.

Study Questions (Answer in writing):

- Who acts as the "gatekeeper" in determining the relevance and reliability of scientific evidence?
- List the criteria used by the Courts to determine if scientific evidence is reliable
- Briefly describe the structure of the Maine court system
- What is a Brady violation?
- Explain the difference between uncertainty of measurement and error rate.
- Define the following terms:
 - Voir dire
 - o Direct examination
 - o Cross examination
 - o Redirect
 - o Chain of custody
 - o Expert witness

Exercises:

- The trainee will prepare an up-to-date curriculum vitae
- The trainee will prepare written answers to the following questions:
 - o Could you please introduce yourself to the jury?
 - o How are you currently employed?
 - o How long have you been employed by the State of Maine as a Chemist in the Tox section?
 - o What are your duties/responsibilities as a Chemist in the Toxicology section?



- O Are you certified by the State of Maine, Dept of Health and Human services (DHHS) as a Blood Alcohol analyst?
- o What exactly does that mean, and how do you gain such certification?
- o Can anyone gain this certification?
- o Did you have to complete any training before getting this certification?
- o Can you tell us about that training? Who provided it and what did it cover?
- o How were you employed before this position?
- o Can you tell us a little about your educational background?
- o Are you a member of any professional societies?
- O You already described your basic duties/responsibilities a Chemist in the Toxicology section, but can you elaborate on the type of testing that is conducted with the toxicology section? Meaning how does the lab receive samples? From whom? How do they get to the lab and how do they move thru the lab?
- o And you, on behalf of the laboratory, issue a report that details your findings?
- O Does anyone in the lab review your findings before they are released? Who and how is that accomplished?
- o Do all blood samples submitted to the lab contain alcohol?
- What happens if you examine a sample that is supposedly from an OUI suspect, but you don't find any alcohol?
- o Can you tell us what a proficiency test is?
- o Are you required to work 'proficiency tests'?
- o How often?
- o Are results from proficiency tests reviewed? By whom?
- o What happens if you 'fail' a proficiency test?
- o Have you ever failed a proficiency test?
- o Can you explain (in lay-terms) what a GC-Headspace system is, and how it determines the amount of alcohol in a blood sample?
- o How do you know the result is just ethanol and not other chemicals and ethanol?
- The trainee will satisfactorily complete a mock trail scenario(s), to include the following sections:
 - Qualifying information
 - o Recognition and introduction of physical evidence
 - O An explanation of the examinations conducted and their results, using visual aids and presentation materials as applicable
 - Conclusions and interpretations



Checklist: Criminal and Civil Law, Procedures, and Testimony

1. The trainee has observed at least 2 courtroom testimonies of a court qualified examiner within HETL, and at least 1 being in the toxicology discipline. Trainee has discussed the testimony with the analyst who testified after the case is completed.

	Date	DA	witness	Detendal	nt	Court/Cnarge
	Trainee		Training Coordinato	r Dat	e	
2.	The trainee	has satisfactorily	prepared an up-to-dat	e curriculum v	vitae.	
	Trainee	·	Training Coordinato	r Dat	e	
3.	The trainee	has prepared satis	sfactory written answe	rs to the outli	ned qu	nalifying questions.
	Trainee		Training Coordinato	r Dat	e	
4.	The trained questions.	e has provided s	atisfactory written/or	al (circle one	e) ans	wers to the study
	Trainee		Training Coordinato	r Dat	e	
5.			Forensic Lab Director t malfeasance in gener		lanage	er to discuss Brady
	Trainee		Training Coordinato	r Dat	e	

Blood Alcohol Training Manual: Doc # = 010 Originally issued 5-8-2017 Approved by: Forensic Lab Director – Lauren Niskach Date Revised: 07-19-2019

Page **19** of **20**



`rai	ning Coordinator	Forensic Lab	Director	Date	
AC	Training Manual Completed:	:			
	Trainee	Training Coo	ordinator	Date	
	Criminal and Civil Law, P reviewed.	rocedures, and	d Testimon	y documen	ted, completed, a
	Trainee	Training Coor	dinator	Date	
	<u> </u>				1
	Bullcoming vs. New Mexico				
	Mitchell V. Wisconsin	ij on orgino			
	U.S. Justice Department Police	ev on Giglio			
	Giglio vs. United States	10118			
	Brady Disclosure Practical Guide to Brady Mot	ions			
	Maine OUI Laws				
	Current Rules of Evidence in	Maine			
	Maine Court Systems				
			Complet	ion Date	Trainee Initials
7.	The trainee has completed rec	quired reading/v	viewing of th	e following	materials:
	Trainee	Training Coor	dinator	Date	
	Conclusions and interpretation				
	Examination/results/presentat				
	Qualifications Recognition/ introduction of 6	nyidanaa			
	Onelifications		Complet	ion Date	Trainee Initial

Blood Alcohol Training Manual: Doc # = 010 Originally issued 5-8-2017 Approved by: Forensic Lab Director – Lauren Niskach Date Revised: 07-19-2019

Page **20** of **20**