

Department of Health and Human Services Health and Environmental Testing Laboratory 221 State Street # 12 State House Station

Augusta, Maine 04333-0012 Tel: (207) 287-2727; Fax: (207) 287-6832

TTY: 1-800-606-0215

Detection of E.coli 0157:H7 /NON 0157 Shiga Toxin Producing E. coli (STEC) by Culture, Real Time PCR and Molecular Serotyping

Laboratory Submission Information Sheet

Laboratory Submission Information Sheet	
Reporting of suspect case to Maine CDC	•Yes: Shiga Toxin Producing E. coli (STEC) is a Notifiable Condition. If you have not already done so, please report suspect or confirmed cases to the Maine CDC via the disease reporting line: 1-800-821-5821 (24hrs/day 7 days/week.)
Required To Submit Laboratory Specimen	•YES
Required Information	•Information on requisition must include: suspected organism if known or required, patient name, DOB, date of collection, specimen source or type, or isolate, submitter name and contact information. •Important: all specimens must be labeled with patient name and be accompanied by a HETL requisition.
Specimen Requirements	• Stool samples in Cary-Blair or MacConkey broth • Raw stool • Isolates
Collection Instructions	• Stool: If possible, collect before start of treatment. Specimen may be raw stool, or stool in Cary-Blair or MacConkey broth. Samples in Cary-Blair should arrive at the laboratory within 48 hours of collection. Raw stool samples must be transported to lab within 2 hours of collection. Isolates should be sent on appropriate biological growth media. Ship samples at ambient temperature, raw stool must be shipped cold
Turn Around Time	•RT-PCR results are generally available within 18-36 hrs of specimen receipt. •Molecular serotyping for STEC is 10 business days •Samples requiring culture and isolation may extend turn-around-time by 72 hours.
Unacceptable Conditions	 Specimens received without a form and clinical details or with insufficient material to analyze. Specimens in poorly labeled/unlabeled or leaking containers. Raw specimens not refrigerated or frozen, or not received in the lab within 2 hours of collection
Results	• Shiga toxin NOT detected, Shiga toxin STX 1 and/or STX2 detected, Shiga toxin STX1 and/or STX2 and E. coli O157:H7 detected, Inconclusive for Shiga toxin • STEC serotyping results • Results will only be reported to the submitter listed on the requisition.
Laboratory Testing Fee	•Fee Waived.
Applicable CPT Codes	•87798 (PCR) •87046 (culture-clinical specimens)
Additional Information	•Bacteriology at 207-287-1704 •Virology, Respiratory, Arbovirus, Serology, and STI at 207-287-1722 •Molecular Biology, Sequencing, and Foodborne at 207-287-5769 •Mycobacterium and Rabies at 207-287-1706