



Forensic Chemistry Section

Fiber Analysis Method

1. **Scope**

This document outlines the general analysis scheme for fiber examinations.

2. **Safety**

- 2.1 The examiner will wear a disposable laboratory coat and disposable gloves as needed when handling evidence for fiber examinations.
- 2.2 The examiner will wear cryogenic gloves and protective goggles when handling liquid nitrogen.

3. **Examination**

- 3.1 Fiber examination will generally proceed in the following order until a point where the examiner concludes that questioned and known fibers are dissimilar. If no such conclusion is reached, the examiner will determine the fibers to be similar.
 - macroscopic examination
 - microscopy (stereomicroscope, compound microscope, polarized light microscope, comparison microscope)
 - microspectrophotometry
 - FTIR microspectroscopy
- 3.2 The analyst will note the macroscopic qualities, which may include but are not limited to amount of sample, texture, color, etc.
- 3.3 A microscopic examination is performed using a microscope with a minimum magnification of 100X and a Fiber Worksheet will be completed.
- 3.4 Microscopic characteristics of natural fibers to be considered will be type of fiber (animal, vegetable, etc.), color, color distribution, and thickness.
- 3.5 Microscopic characteristics of synthetic fibers to be considered will be color, color distribution, thickness, luster and degree of luster, and optical cross section. A polarized light microscope will be used to determine optical properties including parallel and perpendicular relative refractive indices, sign of elongation and birefringence.
- 3.6 The fibers will be examined side by side on a comparison microscope. The examiner will consider color, thickness, cross section, and overall appearance of each fiber.
- 3.7 The comparison will be documented in the examiner's notes. The examiner may choose to also include color printouts or photographs.



Forensic Chemistry Section

Fiber Analysis Method

4. Cross Sectioning

- 4.1 If enough sample is present, and the examiner prefers, a physical cross section may be made of synthetic fibers.
- 4.2 The examiner will place a small (approximately 1cm²) piece of polypropylene on a glass slide.
- 4.3 A fiber will be placed on the polypropylene and covered with another piece of polypropylene of comparable size.
- 4.4 The polypropylene will be covered with another glass slide.
- 4.5 The slides will be heated enough to melt the polypropylene, but being careful not to melt the fiber.
- 4.6 One of the glass slides will be removed.
- 4.7 The polypropylene will be loosened from the other slide.
- 4.8 A scalpel or razor blade should be used to cross section the fiber into very thin slivers. A stereomicroscope is useful for this step.
- 4.9 The cross section portions should be mounted for further microscopic examination.

5. Instrumental Analysis

- 5.1 If the fibers have sufficient color, an instrumental comparison of the color of the known and unknown samples will be performed using a microspectrophotometer. Spectra will be included in the case notes.
- 5.2 If sufficient sample is present, synthetic fiber samples will be analyzed for chemical composition using an FT-IR. The fiber may need to be de-mounted as follows:
 - 5.2.1 Under the stereomicroscope, locate the fiber(s) to be removed.
 - 5.2.2 Use a carbide tipped scribe and press a hole in the cover slip just to the side of the fiber.
 - 5.2.3 Add a small amount of solvent (xylene or xylene substitute) to the hole in the cover slip.



Forensic Chemistry Section

Fiber Analysis Method

- 5.2.4 Use forceps to gently chip away the cover slip to reveal the fiber.
- 5.2.5 Briefly wash the fiber in the solvent and allow to dry.
- 5.3 Samples that need interpretation, i.e. subgeneric class determination, will be compared to spectral libraries present in the software.
- 5.4 When a sample is determined to be of a particular fiber class, the examiner may choose to prepare a known fiber sample and obtain a spectrum for comparison purposes. These types of samples will be treated as the other types of known and unknown samples.
- 5.5 Spectra will be included as part of the examiner's case notes including reference spectra used to identify fiber class.