



Forensic Chemistry Section

Fire Debris Liquid Sample Handling Method

1. Scope

This document provides the procedure for analyzing liquid samples submitted for fire debris analysis.

2. Safety

- 2.1 Fire debris will be processed in the laboratory fume hood. The sash will be lowered as low as practicable when handling the evidence.
- 2.2 The examiner will not sniff or smell the evidence.
- 2.3 Disposable laboratory coats and 6 mil or thicker nitrile gloves will be worn when processing the evidence.
- 2.4 When performing flammability testing, the examiner will take precautions to ensure that the flame is contained to only the sample in question.

3. Method

3.1 Examination

- 3.1.1 Note the condition of the container, including markings and seal.
- 3.1.2 Note the number of liquid layers and a description of each layer.
- 3.1.3 Conduct miscibility testing of each layer.
 - 3.1.3.1 Place an aliquot of the sample liquid in clear vial or test tube.
 - 3.1.3.2 Add water.
 - 3.1.3.3 Agitate the vial or tube, then allow liquid to settle.
 - 3.1.3.4 Note if the liquid is miscible. A single layer indicates the liquid is miscible in water. Two layers indicate the liquid is not miscible in water.
- 3.1.4 Conduct flammability testing of each layer.
 - 3.1.4.1 Place a small amount of liquid on end of cotton swab.
 - 3.1.4.2 Apply flame to sample.



Forensic Chemistry Section

Fire Debris Liquid Sample Handling Method

- 3.1.4.3 Note whether the sample is flammable (e.g. ignites readily, minimal char, large flame) or non-flammable (e.g. does not ignite, does ignite but will not sustain flame, significant char).
- 3.1.5 If the liquid is not ignitable, the examiner may choose to prepare the sample for simple headspace extraction and / or passive diffusion headspace extraction.
- 3.1.6 If the liquid is ignitable and suspected of being a light product or an alcohol, a simple headspace extraction will be performed.
- 3.1.7 If the liquid is ignitable and not suspected of being a light product or alcohol, the sample will be diluted to 20 uL/ml of an appropriate solvent. The solvent will be noted in the case notes.
 - 3.1.7.1 Examiner will save a neat sample, the dilute sample, and the processed sample. A carbon strip will be placed in the dilute sample for long-term storage.