



Forensic Chemistry Section

Method for Preparing Ortho-Tolidine Spray

1. Scope

This document outlines the method for preparing ortho-tolidine spray.

2. Safety

- 2.1 Ortho-tolidine is carcinogenic. Precautions must be taken to minimize exposure. Disposable laboratory coats and gloves will be worn while preparing the reagent.
- 2.2 Reagent preparation will occur in a total exhaust fume hood with the sash as low as practicable.

3. Reagent Preparation

3.1 Acetate Buffer will be prepared as follows:

- 25.0 g sodium acetate will be dissolved in 215 ml glacial acetic acid.
- 250 ml reagent water will be added to the solution.
- The above constituent amounts may be varied as long as the final concentrations remain the same (final volume of approximately 50 ml achieved by using 2.5 g sodium acetate, 21.5 ml acetic acid, and 25 ml water).
- The bottle will be labeled, dated, and stored in the refrigerator.

3.2 The ortho-tolidine spray is prepared as follows:

- 3.2.1 0.4 g ortho-tolidine dihydrochloride will be dissolved into 10 ml of acetate buffer in a small erlenmeyer flask using a magnetic stirrer.
- 3.2.2 The solution will be stirred until it appears creamy.
- 3.2.3 0.5 g sodium perborate will be added to the solution and mixed until most of the solids appear dissolved. The solution will be removed from the stirring plate.
- 3.2.4 35 ml collodion will be mixed with 15 ml ethanol/reagent alcohol in a large erlenmeyer flask with a stirring bar on the stirring plate for 3-4 minutes. **DO NOT HEAT.**
- 3.2.5 120 ml ethyl ether will be **SLOWLY** added to the collodion solution with constant stirring to avoid precipitation.



Forensic Chemistry Section

Method for Preparing Ortho-Tolidine Spray

- 3.2.6 The ortho-tolidine/sodium perborate solution will be added to the collodion solution and stirred for another 3-4 minutes. The resulting solution should appear golden.
- 3.2.7 The solution will be tested against a swabbing of known blood. A blue color should begin to develop within 20 - 30 seconds.
- 3.2.8 The solution will be poured into an aerosol spray bottle and labeled with the name of the solution, date made, initials of preparer, and expiration date (one day).
- 3.3 The remaining unused solution is poured into the appropriate waste container.