



# Forensic Chemistry Section

## Paint Analysis Method

### 1. **Scope**

This document outlines the method for analyzing paint evidence.

### 2. **Safety**

- 2.1 Disposable laboratory coats and disposable gloves will be worn as necessary.
- 2.2 The examiner will wear cryogenic gloves and safety goggles when handling liquid nitrogen.

### 3. **Paint Analysis**

- 3.1 Paint examination will generally proceed in the following order until a point where the examiner concludes that questioned and known paint samples are dissimilar. If no such conclusion is reached, the examiner will determine the paint samples to be similar.
  - macroscopic examination
  - stereomicroscopy
  - microspectrophotometry
  - FT-IR microspectroscopy
- 3.2 The examiner will note the macroscopic qualities, which may include but are not limited to amount and form of sample, texture, color, etc.
- 3.3 If the paint sample is on a hard surface such as car door, tool, etc., the examiner will remove the paint sample with a scalpel or razor blade. The examiner will carefully scrape the paint from the item while ensuring the layers remain intact.
- 3.4 The examiner will perform a stereomicroscopic examination. This comparison may include, but is not limited to, color, layer sequence, layer thickness, texture. The examiner will note the results.
- 3.5 An instrumental comparison of the color of significant layers of known and unknown paint samples may be performed using a microspectrophotometer at the examiners discretion. Resulting spectra will be included in the case folder.
- 3.6 The chemical composition of significant layers of questioned and known paint samples will be compared using a FT-IR microspectrometer. Resulting spectra will be included in the case folder.