Colilert® Test Kit

Introduction

Colilert simultaneously detects total coliforms and *E. coli* in water. It is based on IDEXX's patented Defined Substrate Technology® (DST®). When total coliforms metabolize Colilert's nutrient-indicator, ONPG, the sample turns yellow. When *E. coli* metabolize Colilert's nutrient-indicator, MUG, the sample fluoresces. Colilert can simultaneously detect these bacteria at 1 cfu/100 mL within 24 hours even with as many as 2 million heterotrophic bacteria per 100 mL present.

Contents

WP020I contains 20 Snap Packs for 100 mL samples
WP200I contains 200 Snap Packs for 100 mL samples
W050I contains 20 Snap Packs for 50 mL samples
W050BI contains 200 Snap Packs for 50 mL samples

Storage

Store at 2°–30° C away from light.

Presence/Absence (P/A) Procedure

- Add contents of one pack to a 100 mL sample (50 mL for W050I and W050BI) in a sterile, transparent, nonfluorescing vessel.
- Cap vessel and shake.
- Incubate at 35°C ± 0.5°C for 24 hours.
- Read results according to Result Interpretation table below.

Quanti-Tray® Enumeration Procedure (WP020I and WP200I only)

- Add contents of one pack to a 100 mL water sample in a sterile vessel.
- Cap vessel and shake until dissolved.
- Pour sample/reagent mixture into a Quanti-Tray or Quanti-Tray®/2000 and seal in an IDEXX Quanti-Tray® Sealer.
- Place the sealed tray in a 35°C ± 0.5°C incubator for 24 hours.
- Read results according to the Result Interpretation table below. Count the number of positive wells and refer to the MPN table provided with the trays to obtain a Most Probable Number.

Result Interpretation

Appearance	Result
Less yellow than the comparator*	Negative for total coliforms and <i>E. coli</i>
Yellow equal to or greater than the comparator*	Positive for total coliforms
Yellow and fluorescence equal to or greater than the comparator*	Positive for <i>E. coli</i>

Colilert test kit components: WP020I, WP200I, W050I, and W050BI.

- Look for fluorescence with a 6-watt, 365-nm UV light within 5 inches of the sample in a dark environment. Face light away from your eyes and towards the sample.
- Colilert results are definitive at 24–28 hours. In addition, positives for both total coliforms and *E. coli* observed before 24 hours and negatives observed after 28 hours are also valid.

Procedural Notes

- This insert may not reflect your local regulations. For compliance testing, be sure to follow appropriate regulatory procedures.
- Colilert can be run in any multiple tube format. *Standard Methods for the Examination of Water and Wastewater*** MPN tables should be used to find Most Probable Numbers (MPNs).
- If a water sample has some background color, compare inoculated Colilert sample to a control blank of the same water sample.
- If sample dilutions are made, multiply the MPN value by the dilution factor to obtain the proper quantitative result.
- Use only sterile, nonbuffered, oxidant-free water for dilutions.
- Colilert is a primary water test. Colilert performance characteristics do not apply to samples altered by any pre-enrichment or concentration.
- In samples with excessive chlorine, a blue flash may be seen when adding Colilert. If this is seen, consider sample invalid and discontinue testing.
- Aseptic technique should always be followed when using Colilert. Dispose of in accordance with Good Laboratory Practices.

Quality Control Procedures

The following quality control procedure is recommended for each lot of Colilert:

- Inoculate 3 sterile vessels filled with 100 mL sterile water with the following:
 - one with Quanti-Cult™*** *E. coli* or a sterile loop of ATCC**** 25922 or 11775 (*E. coli*)
 - one with Quanti-Cult *Klebsiella pneumoniae* or a sterile loop of ATCC 31488 (total coliform)
 - one with Quanti-Cult *Pseudomonas aeruginosa* or a sterile loop of ATCC 10145 or 27853 (noncoliform)
- Follow the P/A Procedure or Quanti-Tray Enumeration Procedure above.
- Results should match the Result Interpretation table above.

Colilert test kit components: WP020I, WP200I, W050I, and W050BI.

^[*] IDEXX P/A Comparator, catalog #WP104; Quanti-Tray Comparator #WQTC, or Quanti-Tray/2000 Comparator #WQ2KC
^[**] Eaton AD, Clesceri LS, Greenberg AE. Standard Methods for the Examination of Water and Wastewater. American Public Health Association, 1998, Washington, DC.
^[***] Quanti-Cult cultures – IDEXX catalog # WKIT-1001
^[****] American Type Culture Collection 1-800-638-6597

Colilert, Defined Substrate Technology, DST and Quanti-tray are either trademarks or registered trademarks of IDEXX Laboratories, Inc. in the United States and/or other countries. Quanti-Cult is a trademark of Remel Inc. © 2007 IDEXX Laboratories, Inc. All rights reserved. Manufactured under one or more of the following U.S. patents: 4,925,789; 5,429,933; 5,518,892; 5,610,029; 5,620, 865; 5,620,895; 5,753,456 and 5,780,259. Other U.S. and/or foreign patents issued or pending.

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IDEXX Laboratories, Inc., One IDEXX Drive, Westbrook, Maine 04092 USA