

•No

Reporting of suspect case to Maine CDC •No **Required To Submit** Laboratory Specimen •Information on requisition must include: test requested, patient name, DOB, date of collection, specimen source or type, submitter name and contact information. **Required Information** •Important: all specimens must be labeled with patient name, DOB or other unique identifier, and be accompanied by a HETL requisition. • Preferred: Nasopharyngeal swab, Throat swab, Nasal wash Specimen •Other Acceptable specimens: Lower respiratory specimens (bronchial wash), Nasal aspirate, Sputum Requirements •ONLY synthetic-tipped swabs (Polyester, Rayon, Dacron) with plastic shafts should be used. • Specimens should be submitted in viral transport media provided in the kit. •Store specimens at 2-8°C until ready for shipping. Specimens should be submitted to lab as soon as **Collection Instructions** possible. •Use wet ice or cold packs to maintain 2-8°C during delivery. • PCR results are generally available within 1 - 3 business days of specimen receipt at HETL. Turn Around Time • Specimens received without a form and clinical details or with insufficient material to analyze. Unacceptable • Specimens in poorly labeled or leaking containers or unacceptable specimen source will not be tested. Conditions • Do not use wooden shaft or calcium alginate swabs Do not freeze specimens for PCR testing •Adenovirus DNA Detected or Adenovirus DNA Not Detected • Results will only be reported to the submitter listed on the requisition . Results •\$110 Laboratory Testing Fee •87798 Applicable CPT Codes •Bacteriology at 207-287-1704 • Virology, Respiratory, Arbovirus, Serology, and STI at 207-287-1722 Additional Information

• Molecular Biology, Sequencing, and Foodborne at 207-287-5769

• Mycobacterium and Rabies at 207-287-1706

Detection of Adenovirus by Real Time RT PCR

Laboratory Submission Information Sheet