



Department of Health and Human Services
Health and Environmental Testing Laboratory
221 State Street
12 State House Station
Augusta, Maine 04333-0012
Tel: (207) 287-2727; Fax: (207) 287-6832
TTY: 1-800-606-0215

Sudan ebolavirus Viral Hemorrhagic Fever *Laboratory Submission Information Sheet*

Reporting of suspect case to Maine CDC	<ul style="list-style-type: none">• Yes. This is a Notifiable Condition. If you have not already done so, please report suspect or confirmed cases to the Maine CDC via the disease reporting line: 1-800-821-5821 (24hrs/day 7 days/week.)
Required To Submit Laboratory Specimen	<ul style="list-style-type: none">• Influenza Real-Time PCR, Malaria PCR, and Blood Culture (two sets) may be ordered once.• Ebola Real-Time PCR must have Federal CDC approval before ordering. Medical epidemiology at Maine CDC will coordinate with the attending physician on when to order. Patients must be symptomatic for 72 hours before specimen is collected.
Required Information	<ul style="list-style-type: none">• Information on requisition must include: test requested, patient name, DOB, date of collection, specimen source or type, submitter name and contact information.• Important: all specimens must be labeled with patient name and be accompanied by a HETL requisition.
Specimen Requirements	<ul style="list-style-type: none">• see page 2
Collection Instructions	<ul style="list-style-type: none">• see page 2
Turn Around Time	<ul style="list-style-type: none">• see page 2
Unacceptable Conditions	<ul style="list-style-type: none">• Specimens received without a form and clinical details or with insufficient material to analyze.• Specimens in poorly labeled or leaking containers. Do not use wooden shaft or calcium alginate swabs for Influenza PCR. Samples stored at inappropriate temperatures.
Results	<ul style="list-style-type: none">• Results reported will depend upon the testing performed• All results will be reported only to submitter as stated on requisition via mail or fax.
Laboratory Testing Fee	<ul style="list-style-type: none">• All fees associated with testing will be waived by HETL as long as contact has been made with federal and Maine CDC for prior approval.• \$0
Applicable CPT Codes	<ul style="list-style-type: none">• Influenza - 87502, 87503• Ebola Virus PCR - 87798• Bacterial Culture - 87040
Additional Information	<ul style="list-style-type: none">• HETL Bacteriology Lab at 1-207-287-1704. HETL Virology Lab at 1-207-287-1722. HETL LRN-B Program Advisor at 1-207-287-5679. Disease reporting at 1-800-821-5821



Department of Health and Human Services
 Health and Environmental Testing Laboratory
 221 State Street
 # 12 State House Station
 Augusta, Maine 04333-0012
 Tel: (207) 287-2727; Fax: (207) 287-6832
 TTY: 1-800-606-0215

Test	Specimen Source	Volume	Specimen Container	Storage Conditions	Sample Stability	Turn Around Time	Shipping
Influenza Real-Time PCR	Preferred: Nasopharyngeal swab/Throat swab in viral transport media, Nasal wash Other Acceptable specimens: Lower respiratory	Minimum 1.0ml synthetic-tipped swabs with plastic shaft only	viral transport media	2 – 8 °C	6 hrs	< 8hr	IATA Category A
Malaria PCR	Whole Blood	Minimum 0.5 ml	EDTA (purple top) tubes	2 – 8 °C	6 hrs	< 8hr	IATA Category A
Blood Culture	Whole Blood	See mfg protocol	Facility Specific Blood Culture Bottle	Room Temperature	6 hrs	Variable, Maximum of 7 days	IATA Category A
*Sudan ebolavirus Biofire Film Array NGDS Warrior Panel	Whole Blood	Adult: minimum of 4 mL whole blood Pediatric: minimum of 1 mL whole blood	EDTA (purple top) tubes	Room Temperature	24 hours	Mass PHL (Variable)	IATA Category A
				2 – 8 °C	7 days	Mass PHL (Variable)	IATA Category A
				Dry ice	n/a	Mass PHL (Variable)	IATA Category A

***Collecting Specimens for Ebola Real-Time PCR Testing**

1. Specimens should be obtained when a patient meets the criteria for person under investigation (PUI) including patients with clinical signs, symptoms, and epidemiologic risk factors for Ebola virus disease. If the first specimen is obtained 1-3 days after the onset of symptoms and tests negative and the patient remains symptomatic without another diagnosis, a later specimen is needed to rule-out Ebola virus infection.
2. Staff who collect specimens from PUIs should wear appropriate PPE and should refer to Guidance on Personal Protective Equipment To Be Used by Healthcare Workers During Management of Patients with Ebola Virus in U.S. Hospitals, Including Procedures for Putting On (Donning) and Removing (Doffing).
3. Blood must be collected in plastic collection tubes. Do not transport or ship specimens in glass containers or in heparinized tubes.
4. Do not separate and remove serum or plasma from the primary collection container.
5. Samples may either be shipped from HETL to the Massachusetts Public Health Laboratory, Wadsworth Public Health Laboratory, or Federal CDC.