



Department of Health and Human Services
Maine Center for Disease Control and Prevention
286 Water Street
11 State House Station
Augusta, Maine 04333-0011
Tel.: (207) 287-8016; Fax: (207) 287-9058
TTY Users: Dial 711 (Maine Relay)

Tel. (207) 287-2070

Drinking Water Program

Fax (207) 287-4172

INSTRUCTIONS FOR COLLECTING COMPLIANCE WATER SAMPLES FOR Inorganics

1. Collect sample at the sample point preprinted on the sample form.
2. Remove any screen, hoses and aerators from end of faucet.
3. Run COLD water for 5 to 10 minutes to make sure water has not been sitting for a long time in pipes or tanks and is fresh from the well.
4. Turn water down so it does not splash.
5. Uncap bottle and fill to shoulder with water. Screw the cap on tightly.
6. Fill out the required information on sample form. **SIGN THE FORM!** The sample will not be processed unless the form has a signature!
7. Return bottle with completed paperwork to lab.
SAMPLE MUST ARRIVE AT LAB WITHIN 48 HOURS OF COLLECTION!

**IF YOU HAVE ANY QUESTIONS....PLEASE CALL THE MAINE DRINKING WATER PROGRAM
AT (207) 287-2070 OR YOUR LABORATORY**